



# SwiftDx

## Mycoplasma Detection Kit


### NEW PRODUCT LAUNCH! Order today at [swiftDX.co.uk](http://swiftDX.co.uk)



## Stop mycoplasma contamination in its tracks!

No sample processing or centrifuging of potentially contaminated samples

### Simple & quick test procedure



1 Add 2µl of cell culture media directly to our reaction mix

2 Run reaction in thermal cycler

3 Run reaction on test strip

## Keep possible contamination contained

Quick & simple mycoplasma detection to stop further contamination of your precious samples

Specific & Reliable

Quick Results

No Sample Processing

### Mycoplasma contamination

Mycoplasma contamination of samples can invalidate results, costing drug discovery companies and academic research organisations time, money and reputation. Given the problems it can cause, life science companies are adopting zero-tolerance policies for mycoplasma.

*SwiftDx Mycoplasma Detection Kit allows rapid, routine testing so companies can stop mycoplasma contamination quickly in-house at source.*



### Contact us

[orders@swiftDX.co.uk](mailto:orders@swiftDX.co.uk)

[SwiftDx.co.uk](http://SwiftDx.co.uk)

25 Cambridge Science Park  
Cambridge  
CB4 0FW

